

Safety and Performance of MR-Conditional Pacing Systems with Automated MRI Mode at 1.5 and 3 Tesla

RESULTS FROM THE CAPRI STUDY

Dacher, JN., Langguth, P., Adam, D. et al. Safety and performance of MR-conditional pacing systems with automated MRI mode at 1.5 and 3 Tesla. *Eur Radiol* (2023). <https://doi.org/10.1007/s00330-023-09650-9>


Background & objective

28.5% of Pacemaker patients will have at least 1 MRI scan within 4 years post-implantation.

→ Clinical need for development of MRI conditional pacing systems.

Objective: Evaluate at 1.5T and 3T MRI the safety and performance of trademarked ENO™, TEO™ or OTO™ pacing systems with automated MRI Mode and the image quality of non-enhanced MR examinations.

Methodology



➢ **CAPRI is an interventional, prospective, multicenter,** open label, two arms parallel non-comparative non-randomized study conducted in **Europe and Asia-Pacific** to confirm the clinical safety and performance of the ENO™, TEO™ or OTO™ MRI conditional single or dual chamber rate response pacemakers with VEGA™ pacing leads in the MRI environment.

AUTOMRI™

➢ MRI mode of pacemaker was activated and set to 'Auto' by the cardiologist, allowing **automatic switch to programmed asynchronous pacing mode when the magnetic field of the MRI scanner is detected** and **back to initial configuration, 5 minutes after the patient gets away from magnetic field.** Patients eligible at the inclusion visit underwent an MRI scan at either 1.5 or 3 T.

Results

MicroPort MR Conditional Pacemakers ENO™, TEO™ and OTO™ are SAFE under 1.5 and 3 Tesla MRI scan.

- **Freedom from MRI-related complications** at 1-month post-MRI was **100% in both 1.5 and 3T arms** (both $p < 0.0001$).
- The stability of pacing capture threshold was respectively at 1.5 and 3T [ventricular: both 100% ($p < 0.0001$); atrial: 98.9% ($p = 0.001$) and 100% ($p < 0.0001$)].
- The stability of sensing was respectively at 1.5 and 3T [ventricular: 100% ($p < 0.0001$) and 99.1% ($p = 0.0001$); atrial: 100% ($p = 0.0001$) and 96.9% ($p = 0.01$)].

100% AUTOMRI™ mode proper functioning.

- **All devices switched automatically** to the programmed asynchronous mode in the MRI environment and resumed to normal mode, post MRI exam.
- Even if artefacts were noticed in a subset of examinations, overall **interpretability was preserved.**

Conclusion

This Study demonstrated the **clinical safety and performance of the ENO™, TEO™ or OTO™ MRI conditional pacing systems** at 1-month post-MRI with no scan exclusion zone under **1.5 and 3T**. This study confirmed **the proper functioning of the automated MRI mode.** 100% of MRI images were assessed as interpretable and image quality was assessed as good or excellent in > 90% of the scans.